

DATASHEET

The Pediatric ALS Infant Complete features all important PALS procedures like 4-lead ECG monitoring, training of i.v. and i.o. techniques as well as realistic CPR and airway management training.

The cranial sutures, fontanels and the scalp vein for palpation enable a realistic training.

The intravenous access can be trained on arm and hand and both legs let you train a realistic intraosseous needle insertion with aspiration of bone marrow. The skins and bones are replaceable.

The manikin allows airway management training inclusive Jaw thrust and Sellick's maneuver. It accepts ET tubes up to a size of 3.5 mm and can be intubated orally and nasally. In addition, it accepts NG tubes for feeding procedures.

ECG Monitoring

With the Interactive ECG Simulator 17 pediatric heart rhythms can be recreated, as well as simulated pacing or cardioversion.

It can be used with manual, semi-automatic or automatic defibrillators.

Pediatric ALS Trainer Complete with Interactive ECG Simulator



KEY BENEFITS

- Oral and nasal intubation
- NG tube
- Airway management & CPR training
- Jaw thrust and Sellick's maneuver
- I.V. access, hand and arm
- I.O. access, legs bilateral
- 4-lead ECG monitoring
- Brachial pulse, manual
- Cranial sutures & fontanels
- Scalp vein for palpation

ORDERING SPECIFICATIONS

Item no.	
101-091U	Pediatric ALS Trainer Complete with Interactive Simulator Fullbody manikin, Interactive ECG Simulator, i.v. arm/hand, i.o. legs, pulse bulb, blood powder, soft carry bag
101-092U	I/O Bones & Blood Refill 12-Pack
101-093U	Leg Skins 4-Pack (Pediatric ALS Trainer)
101-094U	Arm/hand Skin 1-Pack (Pediatric ALS Trainer)
101-095U	Arm Vein Replacement Kit

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com