

# SINGLE-USE, NEW POSSIBILITIES



Ambu® aScope™ 4 Cysto

**Ambu**


# aScope™ 4 Cysto OFFERS CONSISTENT QUALITY

Do you notice the gradual deterioration of the image and bending quality of your reusable flexible cystoscopes?

aScope 4 Cysto offers clear and sharp imaging, smooth manoeuvrability and consistent quality. There is no deterioration affecting image or bending quality, and you don't have to worry about damaging your expensive cystoscope. Because with aScope 4 Cysto, you get a brand-new cystoscope for every procedure.

High bending angles of 210°/120° and the varied stiffness of the insertion cord help you manoeuvre and navigate smoothly in the urethra and bladder.

You can advance and completely retroflex the cystoscope to inspect the bladder neck - with or without a biopsy forceps.



## A STRAIGHTFORWARD PROCEDURE

The aView 2 Advance monitor is simple to place in a comfortable position close to the patient.

It's easy to take images and videos and store them to go over the results with the patient in peace and quiet after the procedure or to confer with colleagues.

# aScope™ 4 Cysto IS ALWAYS AVAILABLE AND PORTABLE

Have you ever had to delay a procedure because you didn't have a flexible cystoscope ready when you needed it?

## Make your workday smoother

aScope 4 Cysto is always ready when you are and easy to take with you. Accessing an available cystoscope is no longer a barrier to increasing procedure volume and productivity. Delays that cause your patients to feel disappointed or anxious become a thing of the past.

From the outpatient clinic to inpatient consult procedures, you can take the lightweight single-use cystoscope and portable monitor with you to the Emergency Room or Intensive Care Unit under your arm.

You don't have to locate and move a heavy tower around. And when you're done, you simply dispose of it; there's no more hassle. The result? It's easier to plan and manage your day and to deal with scheduled procedures and unscheduled in-house consults.

## A LIGHTWEIGHT FLEXIBLE CYSTOSCOPE SOLUTION

The aScope 4 Cysto single-use flexible endoscope is always available when you need it.

It weighs just 158 grams, and the compatible aView 2 Advance HD monitor weighs just 2.7 kilos.



## aScope 4 Cysto flexible single-use cystoscope

Item code	Product	Quantity
601001000 (Standard Deflection)*	aScope™ 4 Cysto	5 pcs
600001000 (Reverse Deflection)**	aScope™ 4 Cysto	5 pcs

## aView 2 Advance monitor

405011000***	aView 2 Advance monitor	1 pcs
--------------	-------------------------	-------

## Accessories

K12-11412****	4-way Stopcock. 50 psi (3.5 bar). Swivel Nut.Large Bore	50 pcs
K12-11413****	Hemostasis Valve.AccessPLUS™. Large Bore	25 pcs

\* Standard deflection, commonly used in the US: Moving the control lever downward will make the tip bend down (posteriorly). Moving it upward will make the tip bend up (anteriorly).

\*\* Reverse deflection, commonly used in the EU: Moving the control lever downward will make the tip bend up (anteriorly). Moving it upward will make the tip bend down (posteriorly).

\*\*\* Item number differs by region. Please refer to local material and/or contact Ambu to obtain the full item number for placing orders.

\*\*\*\* Not available in all countries.

# aScope™ 4 Cysto ELIMINATES REPROCESSING AND COSTLY REPAIRS


Have you considered how your clinic could benefit from an optimised workflow?

aScope 4 Cysto enables you to cut costs, such as cleaning supplies, maintenance and repairs (and the associated long-term service contracts). In fact, it may actually save money for your hospital as a whole because it improves your workflow and reduces the time and costs needed to clean reusable cystoscopes.

Also, because it is 100% sterile, aScope 4 Cysto eliminates the risk of cross-contamination.

You can assure your patients that you will be performing the procedure with a single-use sterile cystoscope opened straight from the package just for them.

For more information, please visit [ambu.com/Cysto](http://ambu.com/Cysto)



## REALLOCATE TIME AND RESOURCES

If you process your cystoscopes locally, you can repurpose reprocessing facilities for extra throughput capacity.

You can also free nurses or staff from the hassle and exposure to chemicals when they clean the cystoscope after the procedure. And you can reallocate their time to important tasks like assisting on more procedures.



# TAKE CONTROL OF YOUR SCHEDULE AND BE MORE PRODUCTIVE

What if there were a way to simplify and optimise your workflow and perform more procedures without compromising on quality?

aScope™ 4 Cysto is a single-use flexible endoscope solution that gives you a way to take control of your schedule and be more productive - without compromising on the quality of your work.

It offers consistent quality because you get a brand-new cystoscope for every procedure. It has the image quality and bending performance you need to perform your cystoscopies confidently. In addition, it is always available and portable making it easier to manage your schedule and deal with in-house consult procedures. Finally, it eliminates the need for reprocessing, simplifies workflow, frees up resources and enables you to treat more patients.

## FOCUS ON SINGLE-USE ENDOSCOPY RESULTS IN EXCELLENCE

In 2009, Ambu launched the world's first single-use flexible bronchoscope, Ambu® aScope™, and in 2020, over one million of our single-use endoscopes were used in more than 6000 hospitals globally. The aScope 4 Cysto builds on this decade of experience and know-how.



#### **A HISTORY OF BREAKTHROUGH IDEAS**

Ambu has been bringing the solutions of the future to life since 1937. Today, millions of patients and healthcare professionals worldwide depend on the efficiency, safety and performance of our single-use endoscopy, anaesthesia, and patient monitoring & diagnostics solutions.

The manifestations of our efforts have ranged from early innovations like the Ambu® Bag™ resuscitator and the Ambu® BlueSensor™ electrodes to our newest landmark solutions like the Ambu® aScope™ - the world's first single-use flexible endoscope. Moreover, we continuously look to the future with a commitment to deliver innovative quality products that have a positive impact on your work.

Headquartered near Copenhagen in Denmark, Ambu employs approximately 4,600 people in Europe, North America and the Asia Pacific.

#### **AMBU TAKES A STAND ON THE ENVIRONMENT**

As a company founded and headquartered in Denmark, Ambu is inspired by the positive environmental example of other Danish companies. Furthermore, as a global supplier of single-use medical devices and the largest and most innovative supplier of single-use endoscopes, we feel a responsibility to forge the way with sustainable practices and targets.

#### **AMBU ENVIRONMENTAL TARGETS**

- 100% phthalate-free products by 2020 - achieved
- 95% PVC-free for products released after 2025 - on target
- All packaging reusable and biodegradable by 2025 - on target

**For more information, please visit [ambu.com/cysto](https://www.ambu.com/cysto)**

Ambu A/S  
Baltorpbakken 13  
2750 Ballerup  
Denmark  
Tel.: +45 7225 2000  
[ambu.com](https://www.ambu.com)